


BIOACT® 105 Liquid / wipes



Features:

- ✓ Replaces IPA, MEK, and acetone, High flash point (105F),
- ✓ Low toxicity Fast drying time,
- ✓ Reduced solvent usage
- ✓ High-purity solvent, leaves substrate with very low surface residue
- ✓ Available: (liquid) in pails or drums
(wipes) box of 6 canisters (70 wipes/canister)

Uses:

- ✓ Multi-purpose soak, hand-cleaning, or pre-saturated wipes solvent cleaner.
- ✓ Safe on: ferrous and non-ferrous alloys such as aluminum, copper and copper alloys, magnesium, nickel, steel, stainless steel, titanium & zinc, compatible with many polymeric substrates and elastomers.
- ✓ Removes: adhesives, asphalt, buffing compounds, corrosion inhibitors, coolants, fingerprints, greases, hard water stains, hydraulic fluids, inks, mold release compounds, paint overspray, particulate matter, polishing compounds, sealants, silicones, tars & waxes.

Aerospace approvals:

Bell Helicopter (No 1915)
CFM (CP# 2719)
Gulfstream Aerospace (GAS 115H1)
IAE (CoMat 01-506/A)
Rolls-Royce (MLC204, OMat 1/257S)

Bombardier (BAPS 180-009)
GE Aviation (C04-251)
Honeywell Aerospace (C5409)
Pratt & Whitney (PMC 8929)
United Airlines